

FlowAnalyser™ Set V

Biomedical test set

Check and calibrate all types of gas flow and pressure-producing devices in a very short time. Create your individual report for each test.

The biomedical test set "FlowAnalyser™ Set V" is a certified, high precise measuring system delivered with all needed component parts.



Order number: 300.116.003



FlowAnalyser™ PF-300

The FlowAnalyser measures flow, pressure, temperature, humidity and O2 concentrations bi-directionally. The device for measurements and calibration combines a simple, intuitive multilingual user interface with the highest precision.

FlowLab™ Software

FlowLab is the ideal software solution which has an impressive simple menu navigation.

Biomedical test set "FlowAnalyser™ Set V" to test and verify:

- ✓ Ventilators CPAP / Bilevel
- ✓ Ventilators ICU
- ✓ Ventilators Infant
- ✓ Ventilators High Frequency
- ✓ Blood pressure analysers
- ✓ Oxygen concentrators
- ✓ Vacuum pumps
- ✓ Spirometers
- ✓ Pipe gases

Biomedical test set

"FlowAnalyser™ Set V" includes

FlowAnalyser™ PF-300 / PF-301 / PF-302, SmartLung™ Adult, FlowLab™ Software, adapterset, bacterial filter, USB cable, power supply unit, carrying case, user manual